



Bullard EVA

The Protection Anesthesiologist Need...Everyday!

EVA represents the latest in cutting-edge Powered Air-Purifying Respirator (PAPR) technology and engineering. With many innovative features, such as Active Air Flow Technology™, single HEPA filter, two speed air flow and more...the Bullard EVA is truly revolutionary!

“Anesthesiologists should...
become familiar with
alternative respiratory protective
devices such as the powered
air-purifying respirator(PAPR).

The PAPR provides a high flow rate
of HEPA-filtered air and offers greater
respiratory and contact protection
than the N95.”

*The American Society
of Anesthesiologists,
July 2009 Issue*



Special Features

VISUAL FUEL GAUGE: The easy to read fuel gauge displays remaining battery capacity at a glance. User can check battery status at any time while using EVA.

ACTIVE FLOW TECHNOLOGY™: Another innovation from Bullard...Active Flow Technology automatically responds to user's need for more or less air flow. This intelligent system is continuously working to maintain constant air flow, regardless of filter type, hood type, filter loading, or battery capacity.

INTEGRATED ALARMS: Two distinct audible alarms notify user of low battery and low air flow. The alarm sound travels through the breathing tube directly to the user's ear.

DURABILITY: EVA is powered by a long-lasting 10,000-hour brushless motor. Unlike other motors that can quickly wear out, EVAs motor is designed for long-term reliable use.

SINGLE HEPA FILTER: EVA offers a single HEPA filter for easy use. While other manufacturers' PAPRs may require up to three combination filters, EVA makes it simple, with an easy-to-install single HEPA filter.

FILTER LOCKING TAB: EVA helps ensure positive filter installation with a single click. The filter locking tab notifies the user that the filter is in place with both a visual indication and audible click.



EVA - Evolutionary Air

The EVA PAPR is a completely self-contained, belt-mounted PAPR with a single HEPA filter and a loose fitting facepiece that offers Anesthesiologist advanced particulate protection at a high efficiency level.

(Loose fitting facepiece needed for the use of a stethoscope).

evolve
emergency
PREPAREDNESS

9 Reasons Why More Anesthesiologist Choose EVA

- 1 EVA Hoods require NO fit testing** and have five times the APF compared to a N95 mask in which fit testing is required. Users can choose from more than 15 hood styles, including loose-fitting facepieces, needed for use with stethoscope. (EVA has an APF of 1,000 when used with double bib hood.)
- 2 EVA is the only PAPR with Active Flow Technology™** to ensure a constant air flow which reacts to the user's need for air and works despite a depleting battery. In addition, EVA has a dual speed airflow. Users can toggle to the desired setting (high speed for more air or low speed for less air).
- 3 EVA is easier to DECON.** EVA was engineered with a chemically resistant, sealed & gasketed housing, and closed cell foam belt. EVA also has a shower cap accessory that enables the unit to undergo DECON shower.
- 4 EVA has two distinct audible alarms** that notifies the users of low battery or low air flow. Alarm sound travels through the breathing tube to the user's ear. Low battery alarm alerts users with at least 20 minutes of battery time remaining. Low airflow alerts user when flow falls below 185 Lpm.
- 5 EVA has a durable brushless 10,000-hour motor** for superior reliability. In addition, EVAs brushless battery is one of the quietest available, while providing significantly higher airflows.
- 6 EVA has a battery that lasts 10 hours** and an easy-to-read visual battery fuel gauge. LEDs illuminate to display remaining capacity in 25% increments.
- 7 EVA has a single HEPA filter** for quick and easy installation and a filter locking tab that offers visual and audible indication of positive filter installation. EVA has lower consumable costs vs. blowers that require up to three combination filters.
- 8 EVA is ergonomically designed,** lightweight, and balanced. EVA comes with a contoured belt made of closed-cell, chemically resistant foam shaped to contour to the back for added comfort.
- 9 EVA comes complete** with everything you need to begin using your PAPR including a custom embroidered water resistant storage bag made of durable canvas.

**Call Us Today at
1-888-315-9007
or visit: evolvecorp.com**

EVA Comes Complete with everything you need to begin using your PAPR!



Loose fitting facepiece to allow for the use of a stethoscope



EVA Series Blower with Comfort Belt to provide maximum lumbar support



Single HEPA Filter



Airflow Indicator



10 Hour Battery with Alarm



Standard Breathing Tube

Battery Charger



Custom Embroidered Water Resistant Canvas Bag and Luggage Tag



Evolve Technologies Corporation is committed to the protection and preservation of the environment through the responsible stewardship of the earth's resources. Read more about [Evolve's Green Policy](#) at www.evolvecorp.com. Evolve Emergency Preparedness is a Division of Evolve Technologies Corporation. Evolve, the Evolve logo and the Evolve Green Program logo are trademarks of Evolve Technologies Corporation, 7 Industrial Way, Salem, NH 03079. All other products/product names are trademarks of their respective holders.

Quote printed with permission and courtesy of the American Society of Anesthesiologists, *Advice on Respiratory Personal Protective Equipment for Care of Influenza Patients*, July 2009, Volume 73, Issue 7. All rights reserved. **Read the entire article at www.evolvecorp.com**



The NIOSH approval for the EVA is 21C-0836 for protection against particulates at a high efficiency level [HE]. EVA is ISO 9001 certified.

Printed on FSC Certified Paper.

evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817



7 Industrial Way, Salem, NH 03079
www.evolvecorp.com

3 WAYS TO ORDER

1. Fax completed order form to 603-894-0074
2. Call 888-315-9007
3. Visit evolvecorp.com to contact Evolve

P.O. #

Order Date

BILL TO ADDRESS	
Contact	<input type="text"/>
Company	<input type="text"/>
Address	<input type="text"/>
City	<input type="text"/>
State	<input type="text"/> Zip <input type="text"/>
Tel	<input type="text"/>
E-mail	<input type="text"/>

SHIP TO ADDRESS (Only if different from billing address)	
Contact	<input type="text"/>
Company	<input type="text"/>
Address	<input type="text"/>
City	<input type="text"/>
State	<input type="text"/> Zip <input type="text"/>
Tel	<input type="text"/>
E-mail	<input type="text"/>

ITEM #	QTY.	ITEM DESCRIPTION	UNIT PRICE	SUBTOTAL
BUL-EVAPAN1	<input type="text"/>	<p>Bullard EVA PAPER Including:</p> <ul style="list-style-type: none"> 1 EVA Series Blower with Comfort Belt 1 Single HEPA Filter 1 Battery Charger 1 10-Hour Batter with Alarm 1 Airflow Indicator 1 Loose Fitting Facepiece 1 Standard Breathing Tube 1 Water Resistant Canvas Bag with Luggage Tag (extra charge for Custom Embroidery) 1 Year Warranty (Blower, Battery, Charger) 	\$875.00	
	<input type="text"/>	<p>Custom Embroidery</p> <p><i>Custom embroidering is available up to 22 characters. Please write one letter per box and include spaces where needed by leaving boxes blank.</i></p> <p><input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/></p> <p><i>Custom embroidering will be located on the right end side of the canvas bag. Evolve logo appears on front of bag. Custom embroidering is in Script only.</i></p>	\$25.00	



Payment Options

If you prefer a member of the Evolve staff will contact you regarding payment options. Please include your telephone number.

Card # _____ Expiration Date _____

Name on Card _____ 3-digit Security Code _____

Sub Total \$:

Shipping Cost: **FREE**

Thank you for your order!