

AirBoss Defense HE/OV/AG Cartridge

For ISO Pod | Filters Harmful Gasses & Particulates

This cartridge is a critical component for use with ISOPOD Isolation Systems. Tested against the gasses to NIOSH 42 CFR Part 84 by independent laboratories and has excellent capability in filtering harmful gasses as well as particulates.



It is the only authorized filter for filtering the “inside air” in the ABD ISOPOD. Filters chemicals and particulates that may be off gassed by a patient in the ISOPOD. Tested and proven effective for protecting first responders and medical personnel involved in the transport and treatment infectious diseases such as Ebola and MERS across the use and around the world. Has a 5 year shelf life when stored in the original poly bag.

Who is it for?

Anyone using an ABD ISOPOD.

Specifications	
Part Number	ABD-S-70022-F01
Shelf Life	5 Years (polybag)
Diameter	4.57 inches (116 mm)
Overall Height	3.50 inches (90mm)
Approx. Weight	12.5 oz (355g)
Connection EN	148-1 (40mm) NATO thread
Body Material	Noryl® GFN1
Pressure Drop	38 mm H ₂ O (nominal) measured at airflow rate of 85 lpm
Protections	Particulates – High Efficiency Filter Ammonia – AM Chlorine Dioxide – (CD) Chlorine – (CL) Formaldehyde – (FM) Hydrogen Chloride – (HC) Hydrogen Fluoride – (HF) Hydrogen Sulfide – (HS) Methylamine – (MA) Organic Vapors – (OV) Sulfur Dioxide – (SD)