

Evolve Infectious Disease Body Bag System

3-Layer System, Meets CDC Guidelines for Safe Handling of Human Remains

The Evolve Infectious Disease Body Bag Systems is CDC compliant and protects against the postmortem spread of infectious disease at the site of death, prior to transport, during transport, at the mortuary, and during final disposition of remains.

Each 3-Layer Infectious Disease Kit Contains:

Outer Bag

Size: Extra Large (XL) Adult Size \ Length: 98" \ Width: 48"

Material: Thickness 500 micron (20 mil) \ White PEVA (polyethylene vinyl acetate) \ Chlorine-Free Material Suitable for Cremation or Burial \ Impermeable & Non-Porous \ Provides Containment of Body Fluids & Blood Borne Pathogens (zippers are not leak proof).

Construction: Heat Sealed, Leak Proof seams

Zipper: Curved, Black U-shape Zipper and 2 Zipper Pulls with Tie Ribs for Easy Manipulation while Wearing Protective Gloves

Identification: Integrated Clear Label Pouch w/ ID Tags

Handles: Twelve (12) Built-in Padded Handles

Capacity: 800 lbs

Midde Layer

Size: Standard Adult Size \ Length: 95" \ Width: 40"

Material: Thickness 200 micron (7.9 mil) \ Silver/Green Color \ Chlorine-Free non-PVC Material Suitable for Cremation or Burial \ Impermeable & Non-Porous \ Provides Containment of Body Fluids & Blood Borne Pathogens.

Construction: 2 Sheets heat sealed on one 95" side \ folded & labeled with CDC instructions \ No handles \ No zipper.

Heat Sealer: Required for sealing of this layer

60ga BON

15# LDPE

30ga Foil

6mil Heavy Green PE

MIL-PRF-13K CLASS 1
CADPAK HD100
Sealing Guidelines
Temperature: 435°F
Pressure: 60PSI
Dwell: 3 Seconds

Inner Bag

Size: Standard Adult Size \ Length: 90" \ Width: 36"

Material: Thickness 200 micron (8 mil) \ White PEVA (polyethylene vinyl acetate) \ Chlorine-Free Material Suitable for Cremation or Burial \ Impermeable & Non-Porous \ Provides Containment of Body Fluids and Blood Borne Pathogens (zippers are not leak proof).

Construction: Heat Sealed, Leak Proof seams \ No handles

Zipper: Curved Black U-shape Zipper and 2 Zipper Pulls with Tie Ribs for Easy Manipulation while Wearing Protective Gloves

Identification: Integrated Clear Label Pouch w/ ID Tags.

Packages

Part# EPCKITSEAL6	6-Pack: Contains Six (6) 3-Layer CDC Infectious Disease Kits & Sealer
Part# EPCKITSEAL3	3-Pack: Contains Three (3) 3-Layer CDC Infectious Disease Kits & Sealer
Part# EPCKITSEAL2	2-Pack: Contains two (2) 3-Layer CDC Infectious Disease Kits & Sealer
Part# EPCKIT1	Single Pack: Contains One (1) 3-Layer CDC Infectious Disease Kit , NO Sealer
Part # EPCSEAL	Sealer: One (1) 12" Heavy Duty Heat Crimper (Sealer ONLY)