



WHAT SHOULD YOU EXPECT FROM YOUR DISPLAY?

THE FULL ADVANTAGE OF MULTI-MODALITY VIEWING

In addition to pixel-perfect grayscale imaging, Coronis Fusion 6MP DL brings you the full advantage of color modalities like 3D reconstruction, color MRI and PET-CT imaging. Equipped with the latest IPS-Pro LCD panel and a unique Diagnostic Luminance backlight technology, Coronis Fusion 6MP DL provides an unprecedented brightness of 800 cd/m² and a superb contrast ratio of 800:1. It means that, from now on, grayscale as well as color studies can be read on a single, multi-modality display system. No more need for workflow-interrupting journeys from your PACS workstation to the modality acquisition workstation to interpret color enhanced images.

**Everything you expect and more...
to experience absolute freedom and comfort**



EXPERIENCE UNIFORMITY

With the unique Uniform Luminance Technology, luminance is spread equally over the entire 30-inch screen. Your diagnostic images are DICOM compliant from center to corner. In addition, the intelligent Ambient Light Compensation system assures that Coronis Fusion 6MP DL remains DICOM P14 compliant under all lighting conditions by adapting its brightness relative to the ambient light level in the room.

- Corner to Corner Luminance Uniformity
- Diagnostic Confidence Under Changing Light Conditions

EXPERIENCE VERSATILITY

The uninterrupted wide-screen desktop gives you the freedom to organize your workspace just the way you want it. In addition, color and contrast are not distorted when imaging is viewed from an angle. IPS-Pro delivers a maximum viewing angle of 178° and a typical contrast ratio of 800:1. IPSPro widens the viewing cone in which contrast levels are high enough to detect even the most subtle lesions.

- True Multi-user Viewing
- Perfect View From Every Angle
- Freedom to Organize workspace

EXPERIENCE TRUE PRODUCTIVITY

For the first time, a 30-inch diagnostic color monitor can be used as two uninterrupted 3MP displays or one wide-screen 6MP display. To optimize productivity, the Coronis Fusion 6MP allows the freedom to read CR, DR, PET-CT, Cath, CT, Colonoscopy and MR images side by side on a single wide screen viewport without a disturbing center bezel.

- Side by Side Viewing
- Maximum Efficiency in 2D and 3D
- Dynamic Images Without Artifacts
- Continuous DICOM Conformance
- Automated Quality Assurance

**Call Us Today at
1-888-315-9007
or visit: evolvetechnology.com**

evolve
medical imaging

7 Industrial Way
Salem, NH 03079-2817