



Decontamination Wheelchair

Cost Effective, Does Not Rust or Sustain Bacterial Growth

The Decontamination Wheelchair has a poly resin plastic frame, thereby eliminating rust and chips. It has 24 inch plastic wheels and hand rims with non-marking tires. The Decontamination Wheelchair has 5 inch total locking front casters and a pull out footrest. Standard seat height is 21 inches. The tension adjustable mesh backrest and sling seat keep the Decontamination Wheelchair in top shape. It comes with a molded backrest with a built-in push handle.

Decontamination Wheelchair Specifications

Between Wheelchair Arms	17.5"
Overall Width	29.5"
Seat Depth	18"
Weight Capacity	250 lbs

Options

Support: Torso Belt, X-Harness, H-Harness, Seat Belt

Calf/Footrest: Grip Tape on Footrest

Current Mesh Colors: Almond, Autumn Fern, Grass Cloth, Green, Grey, Navy Blue, Red, Royal Blue, Tropical Foliage, Wild Orchid Black, Yellow

Additional Options: Handbrakes

Material

Tubing & Fittings: Georgia Gulf #7140 special purpose rigid PVC compound

Mesh: PVC coated polyester mesh (machine washable)

