

Evolve Emergency Modular Hospital System



The AI-77 turns the individual patient bed into an isolation zone in a matter of minutes, giving the healthcare worker protection from airborne infectious disease. Surpasses CDC guidelines for infectious disease isolation requirements.




The Rapid Response Emergency Bed is a cost effective solution that provides immediate additional bed capacity for overflow in the emergency department, a field hospital, or mass casualty incident. Quick assembly, durable and easy to decontaminate.



AIRISO PASS (Portable Air Safety Systems) line of equipment meets all of your air scrubbing or medical isolation needs, as well as odor abatement and general air purification requirements. We reach unmatched levels of performance by employing our patented technology of HEPA filtration and UVC irradiation as a portable system on wheels. The HEPA filter is under 24/7 UVC irradiation, neutralizing all bio-aerosol contamination captured by the filter. All PASS units incorporate high flow HEPA filters rated at 99.99% effective at 0.3 microns.



The DRAKE line of portable ECUs offers all-in-one service for air conditioning, heating and dehumidifying for all field applications. We offer many different mounting configurations to include: cart mounted, trailer mounted and skid mounted. All versions available in 60Hz and 50Hz configurations.



The Evolve Emergency Modular Hospital System is completely customizable. The illustration to the left features the Losberger TAG 5A which fits 12 Rapid Response Emergency Beds providing 45 square feet for each patient.

Choose from a variety of environmental control units, HEPA filtration & isolation systems, and bed solutions. All systems feature Losberger TAG Inflatable Shelters that deploy in less than 10 minutes with 2-3 people.

Call Evolve Emergency Preparedness: 1-888-315-9007

Evolve Emergency Modular Hospital System

The Losberger TAG Inflatable Shelters are deployable in less than 10 minutes. Every TAG shelter can be connected to form larger shelters in any direction (side and gable doors). TAG shelters are made of durable fire and UV resistant PVC coated fabric (military grade 22 oz/sqyd), using HF welding to construct the cover and floor into one enclosure.



Losberger TAG Shelter	Square Feet	Dimensions	Bed Capacity	Isolation Filtration Specifications (Optional)	Environmental Control Unit(s) (Optional)
TAG2A	161	L 7'10" x W 18'4" x H 9'2"	4 RRE Beds	PASS 700	2T ECU (AC and Heat)
TAG3A	312	L 16'11" x W 18'4" x H 9'2"	6 RRE Beds	PASS 700	3.5T ECU (AC and Heat)
TAG4A	452	L 24'9" x W 18'4" x H 9'2"	10 RRE Beds	PASS 700 (Need: 2)	5T ECU (AC and Heat)
TAG5A	603	L 32'8" x W 18'4" x H 9'2"	12 RRE Beds	PASS 700 (Need: 2)	3T ECU (AC and Heat) (Need: 2)
TAG6A	840	L 40'6" x W 20'8" x H 10'2"	16 RRE Beds	PASS 700 (Need: 3)	5T ECU and 3.5T ECU (Need: 1 each)

Complete your field hospital with any combination of these options.

- Rapid response emergency beds (adult, bariatric, pediatric, crib style and veterinary versions available)
- Environmental control units: available for all shelter sizes
- HEPA filtration units: PASS systems for air scrubbing exceed recommended air change levels
- Isolation beds: maximize your bed space with individual HEPA filtration bed systems
- Mobile medical systems: complete clinical packages organized for field deployment: triage, wound care, general medicine, PPE, surgical, OB-GYN, etc.
- Complete lighting packages for Losberger Rapid Inflatable Shelters
- Hard flooring packages available for all shelters
- PPE for clinical and support staff: gloves, disposable bio-protection suits, eye protection, respiratory protection
- Decontamination systems: field deployable, single stage shower and multi-lane, mass decontamination systems available
- Mobile water treatment systems for all field operations

Optional Equipment

- Automatic pressure control system
- Internal partition walls
- Internal liners
- External sunshade
- Protective groundsheet
- Lighting & electrical distribution box
- Heating/AC systems
- Hard flooring
- Side door(s)
- Connecting & entrance modules
- Printed logos



7 Industrial Way
Salem, NH 03079-2817
evolvetechnology.com
1-888-315-9007

