



ICC Rapid Response Emergency Bed

Triage and Emergency Management Preparedness System

Fully Assembled
and Patient Ready
in 30 Seconds



0 Seconds



12 Seconds



23 Seconds



30 Seconds

ICC Rapid Response Emergency Bed provides immediate additional bed capacity for triage during a mass casualty incident or routine Emergency Department surge overflow. The Temps Beds (knock-down gurneys) consist of an easy to assemble, knock-down patient gurney and a color coded, visual indicator that provides critical triage information to medical personnel at a single glance.



Special Features

- Provides immediate additional bed capacity
- Assemble in less than 30 seconds
- Utilizes 5" total lock casters
- 3 Safety belts standard
- Designed to hold X-ray film tray
- Tension adjustable patient support sling
- Mesh sling ideal for "hot zone" wash down
- 5 position elevating backrest
- Narrow profile will not block halls
- Comfortable working height



Field Tested for Reliability

ICC Rapid Response Emergency Bed was designed in conjunction with leading experts in the field of trauma/emergency medical care. It has been field tested in rigorous multiple casualty disaster drills in hospitals in several states.



High Quality Construction

ICC Rapid Response Emergency Bed has a tube structure made of the highest quality polyvinyl chloride resin available. It will not chip, dent, crack, peel, fade, stain or sustain bacterial growth. The patient support sling is vinyl coated polyester mesh and is tension adjustable for maximum patient comfort.



**Call Us Today at
1-888-315-9007
or visit:
evolvetechnology.com**

evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817


ICC Rapid Response Emergency Bed Options




Poly-backed **Bio-substance Isolation Mat** for blood and body fluid containment



Isolation Canopy can be snapped into place within seconds and folds for easy storage



CPR Board slides under the patient providing a solid surface for cardiac procedures



Egg Crate and closed cell foam mattress overlays for longer term use



Each Temps Bed can be converted for use in sand and rough terrain with an optional **All Terrain Chassis**



The **Mobile Transport and Storage Cart**, with dust cover holds 5 Temps Beds in the space of 1 and can be easily moved to any location

ICC Rapid Response Emergency Bed Specifications



Swivel Wheels engineered for easy action and long life, allow the Temps Beds to go where no other gurney can go



When used in conjunction with existing protocols, the **Triage Color-coded Visual Indicator** assists in the rapid assessment and classification of patients. (Available in: Green, Yellow, Red, Black. Other colors available upon request)



- Overall length: 75" (6' 3")
- Overall width: 25"
- Light weight easy one-person setup
- Working height: 31"
- Standard weight capacity: hold up to 300 lbs
- Extra weight capacity bariatric model available: holds up to 450 lbs

**Call Us Today at
1-888-315-9007
or visit: evolvetechnology.com**

evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817