



Med Sled Evacuation Basket and Rack System

Safe Pediatric Evacuation • NICU • PICU • New Born • Infant

Evacuation Basket

- Used for evacuation of high risk and NICU infants requiring safe, secure and hands-free carrying
- Allows for high frequency ventilation and medical attention
- Holds IV bags and life-sustaining equipment
- No assembly required



#MSEB	Med Sled Evac Basket 24"L x 13"W x 7"H
-------	---

Evacuation Basket Rack System

- No assembly required, easily deploys in less than 1 minute
- Baskets simply slide on and lock securely to the Rack
- Single Rack hold six Evacuation Baskets (sold separately)
- Eliminates the need to have infant beds at the alternative care site
- Waist-high for easy accessibility to babies
- Durable forged aluminum holds up to 600 lbs
- Can accommodate four Accessory Poles (IV poles sold separately)
- Decon Capable

#MSEB - RACK	Med Sled Evac Basket Rack System 80"L x 39"W x 39"H (without IV Poles)
--------------	--

evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817
Tel: 1-888-315-9007
evolvetechnology.com

Storage for Evacuation Basket and Rack System



Rack System

Compact
Rack System
Folds to
40"L x 9"D x 8"H



Portable
Rack System
Fits in Carrying
Storage Bag



Evac Basket

Holds Up To 6 Baskets and
1 Rack System (All Sold Separately)



#SDEB6R

Storage Device Holds 6 Med Sled Evac
Baskets with Vests and 1 Eac Basket
Rack System 33"W x 8"D x 40"H

- Easy wall mount storage system - simple to deploy
- Compact construction allows for hanging behind doors or in utility closets - out of the way, yet accessible
- Stores up to 6 Baskets and 1 Rack System
- Durabel canvas protects the baskets from containimation

For More Information Call
Evolve Emergency Preparedness
1-888-315-9007
Or Visit: evolvetechnology.com

evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817