



Med Sled

Tactical Evacuation • Vertical Evacuation • Patient Transport

As a first responder, you realize there is no standard situation for an emergency evacuation. Your environments range from: deep woods, burning houses, school dormitories to architectually unstable hospitals or commercial buildings. Med Sled is an easy to use economical transportation device to help make evacuation safe and fast!



Tactical/Hostile Extrication MS30

- Dramatically reduces the effort of moving wounded. Med Sled enables a one-on-one rescue or evacuation
- Hands-free tether around waist or shoulder allows for carrying equipment or supplies
- Patented design provides stabilization and can be used with a backboard or stokes basket
- Utilizes High Density Polyethylene, a more durable material
- Slides easily for safe and immediate evacuation.
- COMPACT: carrying bag with shoulder strap included.
- LIGHTWEIGHT: only 6lbs.
- DURABLE: made of highly durable plastic material that helps protect person while transporting



Vertical Rescue MS36VLR

Med Sled has significant advantages and benefits over the SKED Sled!

FASTER DEPLOYMENT

- No assembly required
- Integrated hoist system for quick lifting (pre-assembled)

LIGHTER

- Weighs less than 17 lbs
- Fewer components than the SKED Sled enabling immediated deployment

STRONGER

- Utilizes High Density Polyethylene, a more durable material
- 10,000 lb rated perimeter tether secures the Warfighter expediting the evacuation process
- Utilizes 4,000 lb rated cross straps with side release buckles for ease of use
- DECON capable

INTEGRATED HOIST SYSTEM

- Vertical/horizontal lift capability
- Hoist system rated to 10,000lbs
- Helicopter hoist straps included with Med Sled
- No assembly required at time of extraction



Patient Transport MS36

- Simple and easy to use
- Stairwell braking system attached to every Med Sled allows for safe and controlled descents
- No lifting! Uses roll and drag method only
- Enables utilization of all available personnel for a one-on-one evacuation
- Allows staff to transport patients 2-3 times their own weight
- Utilizes High Density Polyethylene, a more durable material
- Patented design provides stabilization and can be used with a backboard or stokes basket , can hold IV bags, oxygen and other small devices
- DECON capable

Call Evolve Emergency Preparedness
1-888-315-9007
Or Visit: evolvetechnology.com

evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817
Tel: 1-888-315-9007
evolvetechnology.com



Med Sled

Tactical Evacuation • Vertical Evacuation • Patient Transport




Tactical/Hostile Extrication MS30



Vertical Rescue MS36VLR



Patient Transport MS36



Each Med Sled Comes With Its Own Carrying Bag With an Over-the-Shoulder Strap



Toddler Insert



Infant Inserts for 0-5 Months

Med Sled Can Help Safely Evacuate All Size Patients!



Med sled Is Easy To Store! (4 Sled Storage Wall Unit)



Bariatric Med Sled

Call Evolve Emergency Preparedness

1-888-315-9007

Or Visit: evolvetechnology.com

evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817

Evolve Emergency Preparedness is a Division of Evolve Technologies Corporation. Evolve, the Evolve logo are trademarks of Evolve Technologies Corporation, 7 Industrial Way, Salem, NH 03079. All other products/product names are trademarks of their respective holders.