



Medical Surge Cot

Economical, No Assembly, Easy-to-decon, Does Not Sustain Bacterial Growth



Medical Surge Cot Specifications

Bed Height - folded for storage	11"
Bed Height - In use	29" floor to patient deck
Overall Length	52"
Folded Length	52"
Overall Width	25"
Folded Width	25"
Weight Capacity	350 lbs standard

Options

Mesh Colors: Almond, Autumn Fern, Grass Cloth, Green, Grey, Navy Blue, Red, Royal Blue, Tropical Foliage, Wild Orchid Black, Yellow

Material


Tubing & Fittings: Georgia Gulf #7140 special purpose rigid PVC compound

Mesh: PVC coated polyester mesh (machine washable)

The Medical Surge Cot provides immediate additional bed capacity for isolation, disaster preparedness, post ER surge, field hospitals, mass casualty incident, evacuation and shelters. Cots are sold individually or in a system of 5 cots with a heavy duty aluminum transport and storage cart that has a water proof cover.

The Medical Surge Cot is constructed of the highest quality polyvinyl chloride resin structural tubing which will not chip, dent, crack, peel, fade, stain or sustain bacterial growth. The patient support sling is a vinyl coated polyester mesh and is tension adjustable for maximum patient comfort. The Medical Surge Cot has been developed in conjunction with leading experts in the field of trauma and emergency medical care and has been field tested in rigorous multiple casualty disaster drills in hospitals in several states.

Medical Surge Cots fold flat for easy transport & compact storage. They are fully assembled and patient ready in less than 30 seconds.



evolve
emergency
PREPAREDNESS

7 Industrial Way
Salem, NH 03079-2817
1-888-315-9007
evolvetechnology.com