



MRIpad

► **Increases Patient Comfort & Stability • Lessens Repeat Sequencing • 100% X-ray Lucent**

The MRIpad is an inflatable mattress pad. It was specifically designed to reduce patient call-backs, and to decrease sedation and patient anxiety by increasing patient stability and comfort during MRI scanning. The MRIpad takes seconds to put in place and can be inflated according to patient weight and comfort needs. The patented engineering renders lateral support and comfort, while decreasing pressure points and filling supine-position voids. The MRIpad can be used with patient's of any size and on most MRI coils.



What Customers Are Saying!

"We have less motion artifact and need fewer repeat studies when we use the MRIpad for patients with back pain. It's a great addition to improving patient care."

Peter Hathaway, MD

"I found it to be a great asset and it provided added comfort for patients who have a hard time lying flat. All of the patients were having an L-spine MRI and were in severe pain. Had it not been for this device, I am uncertain as to whether they would have made it all the way through the exam. They each expressed how much the pad helped and how grateful they were that I had offered it to them."

Paulla Foster, RT (MR)



Increase Throughput and Decrease Patient Call Backs!

The patented design of the MRIpad comforts your patient during the scan; significantly decreasing patient motion, increasing throughput and decreasing patient call backs.

Saves Time and Money!

The MRIpad stabilizes the patient's body relieving pressure and discomfort of hard MRI coils. Repeat scan sequences, due to patient motion, are reduced. Pain management interventions can be decreased saving time and money. The MRIpad is cost effective and can either be reused or disposed - reducing MRSA and other transmitted infections.

Benefits You Can Expect:

- Better Quality Scans
- Decreased Call Backs
- Increased Patient Comfort

evolve
medical imaging

7 Industrial Way, Salem, NH 03079-2817

The MRIPad is made for comfort and is MRI Compatible*

► Increases Patient Comfort & Stability • Lessens Repeat Sequencing • 100% X-ray Lucent

In these side-by-side comparisons, you can see that there is no signal loss and no artifacting of the spine between patients with and without the MRIPad.



Scan with MRIPad.




Scan without MRIPad.



Scan with MRIPad.



Scan without MRIPad.



MRIPad Specifications:

- 100% Latex Free
- Durable PVC construction
- Compatible with most magnets and coils
- 10" x 12" x 1/2" packaged
- Pump sold separately

**Call Us Today at
1-888-315-9007 or visit:
evolvetechnology.com**

* The MRIPad is MRI compatible and is MR safe. It will not significantly affect the quality of the diagnostic information, nor will its operation be significantly affected being in the MR environment.

Evolve Medical Imaging is a Division of Evolve Technologies Corporation. Evolve and the Evolve logo are trademarks of Evolve Technologies Corporation, 7 Industrial Way, Salem, NH 03079. All other products and product names are trademarks of their respective holders.

Printed on FSC Certified Recycled Paper. 

evolve
medical imaging
7 Industrial Way, Salem, NH 03079-2817