



Pediatric Emergency Bed

Cost Effective, Easy to Assemble & Decon, Does Not Sustain Bacterial Growth

Pediatric Emergency Beds fold flat for easy transport & compact storage. They are fully assembled and patient ready in less than 30 seconds.



The Pediatric Emergency Bed has all the features of the standard original Rapid Response Emergency Bed. The Pediatric Emergency Bed is a cost effective solution that provides immediate additional bed capacity for surge overflow in the emergency department, a field hospital, or mass casualty incident. The Pediatric Emergency Bed consists of an easy-to-assemble, knock-down pediatric patient gurney, and a color coded visual indicator that provides critical triage information to medical personnel at a single glance.

Options

Current Mesh Colors:

Almond, Autumn Fern, Grass Cloth, Green, Grey, Navy Blue, Red, Royal Blue, Tropical Foliage, Wild Orchid Black, Yellow

The Pediatric Emergency Bed is constructed of the highest quality polyvinyl chloride resin structural tubing which will not chip, dent, crack, peel, fade, stain or sustain bacterial growth. Therefore, it can be easily decontaminated. The patient support sling is a vinyl coated polyester mesh and is tension adjustable for maximum comfort. The Pediatric Emergency Bed has been developed in conjunction with leading experts in the field of trauma and emergency medical care and has been field tested in rigorous multiple casualty disaster drills in hospitals in several states.

Pediatric Beds are sold individually, or in a system of 5 beds with a heavy duty aluminum transport and storage cart and water proof cover.

Pediatric Emergency Bed Specifications

Bed Height - folded for storage	11"
Bed Height - In use	29" floor to patient deck
Overall Length	52"
Folded Length	52"
Overall Width	25"
Folded Width	25"
Weight Capacity	350 lbs standard

Material

Tubing & Fittings: Georgia Gulf #7140 special purpose rigid PVC compound
Mesh: PVC coated polyester mesh (machine washable)

evolve
emergency
 PREPAREDNESS

7 Industrial Way
 Salem, NH 03079-2817
1-888-315-9007
 evolvetechnology.com