

ScanX Quantum Computed Radiography Dependable, Affordable Full Body & Dental Imaging

ScanX Quantum Introductory Package Includes:

- ScanX Quantum
- **Choose either a desktop computer or laptop**
- PACSmart™ DICOM Software with Image Acquisition and 2 DICOM Viewer licenses
- (2) 14" x 17" Imaging Plates
- (2) 14" x 17" Cassettes
- (2) 14" x 17" Plate Protectors
- (1) 10" x 12" Imaging Plate
- (1) 10" x 12" Cassette
- (1) 10" x 12" Plate Protector
- (8) Size 2 Dental Imaging Plates and Barrier Envelopes
- (1) Size 0 Dental Imaging Plate and Barrier Envelopes
- Dental Plate Transfer Box

ScanX Quantum digital imaging systems are easy to use and versatile, enabling the user to process both full body and dental images. ScanX systems are built on well-known techniques used with film, so procedures are familiar. Technicians can obtain high-quality digital images in just three simple steps. 1.) Expose, 2.) Insert, 3.) View. The ScanX Quantum accepts all phosphor plates in sizes 8"x10", 10"x12" and 14"x17", and all dental sizes zero to six. Upgrading to digital with the ScanX Quantum is a wise investment that will enhance the workflow of your practice and improve your bottom line.



- Compact and portable
- Wall-mountable for increased space savings
- Finely detailed digital images help you to see even the smallest details
- Customizable settings enable optimum image quality for your specific application
- Accepts any brand of phosphor storage plates in standard sizes up to 14" x 17"
- Can process dental images in all sizes
- Easy to operate (staff can independently complete revenue generating dental imaging procedures)
- Patented, built-in eraser with manual or auto erase function (plates are immediately ready for re-use)
- Operates in normal room lighting; darkroom and chemicals are eliminated
- Lightweight extruded aluminum frame
- Built-in rechargeable battery option
(useful features for mobile clinics)

**Call Us Today at
1-888-315-9007 or visit:
evolvetechnology.com**

evolve
medical imaging

7 Industrial Way, Salem, NH 03079-2817