



Surge Crib - 3 Crib System

Fully Assembled, Easy to Decon, Does Not Sustain Bacterial Growth



The Surge Crib provides immediate additional crib capacity for surge overflow, evacuation, field hospitals, shelters or a mass casualty incident. The Surge Crib assembles in 60 seconds. The 3 Crib System utilizes fold, stack and store technology and includes a transport and storage cart and dust cover.

The Surge Crib folds flat for easy transport & compact storage. They are fully assembled and patient ready in less than 30 seconds.

The Surge Crib is constructed of the highest quality polyvinyl chloride resin structural tubing which will not chip, dent, crack, peel, fade, stain or sustain bacterial growth. The patient support sling is a vinyl coated polyester mesh and is tension adjustable for maximum patient comfort. The Surge Crib has been developed in conjunction with leading experts in the field of trauma and emergency medical care and has been field tested in rigorous multiple casualty disaster drills in hospitals in several states.



Options

Current Mesh Colors:
Almond, Autumn Fern,
Grass Cloth, Green, Grey,
Navy Blue, Red, Royal
Blue, Tropical Foliage,
Wild Orchid Black, Yellow

Surge Crib Specifications

Overall Length	45"
Overall Width	31"
Working Height	32"
Patient Area	25" x 40"
Top of Side	28"
Curtain to Mattress	

Material

Tubing & Fittings: Georgia Gulf
#7140 special purpose rigid PVC
compound
Mesh: PVC coated polyester mesh
(machine washable)

evolve
emergency
PREPAREDNESS

7 Industrial Way
Salem, NH 03079-2817
1-888-315-9007
evolvetechnology.com