



Mobile Vaccination Systems by Vericor

Complete, Organized and Appropriate Quantities for Effective Patient Care

Just add staff and vaccine! These complete, mobile, organized vaccination systems contain all of the supplies necessary for a successful mobile clinic. Intuitive packaging increases staff efficiency and reduces vaccination times. Scalable modules enable providers to support 500 to 20,000 person campaigns. So for your next off-site vaccination clinic, leave the mess of boxes back at the office and conduct a professional vaccination campaign with a Response-Friendly™ Vaccination System.

Vaccinator 500 Go-Bag (VS-500GB)

- All-in-one packaging for quick and easy transport.
- Two (2) wheels and a telescopic handle support efficient mobility of supplies.
- Includes a rolling, soft-sided medical case, sharps containers, vinyl gloves, bandages, nurse's drug handbook, safety needles, syringes, and more
- SmartBook™ included



evolve
emergency
PREPAREDNESS

7 Industrial Way, Salem, NH 03079-2817
Tel: 1-888-315-9007
evolvetechnology.com



Vaccinator 2500 Go-System (VS-2500GS)

- Four (4) labeled, organized drawers keep supplies ready for immediate use.
- Front lid converts into a fold-out table for an additional work surface.
- Includes a 4-Drawer Medical Case, biohazard waste containers and liners, vinyl gloves, bandages, nurse's drug handbook, safety needles, syringes, spill kits, and more.
- SmartBook™ included.



Vaccinator 5000 System with Refrigerator (VS-5000R)

- Includes a 42" Workstation with built-in refrigerator, biohazard waste containers and liners, vinyl gloves, bandages, nurse's drug handbook, safety needles, syringes, a certified traceable thermometer, and more.
- Built-in refrigerator supports recommended cold chain management practices.
- Three (3) mobile medical cases enable in-the-field replenishment.
- SmartBook™ included.





Vaccinator 5000 System (VS-5000)

- Same system as the VS5000R with refrigerator, but with a 6-Drawer 42” Workstation instead of the refrigerator model.



EZ Refill Kits

- Replenish depleted vaccination systems.
- Simply “place in the case”.
- SmartBook™ Update includes new critical data of supplies.

Model #	Description	Includes
VS-EZ500	EZ 500 Refill	Supplies to vaccinate 500 people.
VS-EZ2500GS	EZ 2500GS Refill	Supplies to vaccinate 2500 people.
VS-EZ5000R	EZ 5000R Refill	Supplies to vaccinate 5000 people.



MC-24 Workstation



MC-24 Workstation Closed



MC-42R Workstation

Workstations

- Front lid converts into a fold-out table for an additional work surface.
- Drawers are constructed with durable stainless steel facilitating easy decontamination.
- Partitioned drawers and tilt bins increase organization and efficiency.
- Easy lock, 4" wheels enable smooth and secure storage and transport.
- The safe, mobile design has a wide wheel base and low center of gravity.
- Lockable for secure storage.
- Alphanumeric labeling for easy identification and management.
- Compact storage.

Model #	Description	Differentiating Features	L x W x H
MC-24	Workstation, 24", Table	3 Drawers, 6 Tilt Bins	24" x 24" x 48"
MC-42	Workstation, 42", Table	6 Drawers, 12 Tilt Bins	42" x 24" x 48"
MC-42R	Workstation, 42", Refrig/Table	3 Drawers, 12 Tilt Bins, Refrigerator*	42" x 24" x 48"
MC-48	Workstation, 48", Table	6 Drawers, 12 Tilt Bins	48" x 24" x 48"
MC-48R	Workstation, 48", Refrig/Table	3 Drawers, 12 Tilt Bins, Refrigerator*	48" x 24" x 48"

*Built-in, lockable, "stand-alone" refrigerator maintains a consistent temperature for safe vaccine storage (as recommended by the CDC).

Thermometer, Certified Traceable (V-TCT)

- Non-toxic, sealed glycolic liquid probe.
- Programmable alarm.
- Records high and low readings.





Cool Cube™ Vaccine Coolers/Transporters

- No ice or electricity required.
- Easy, quick assembly with a single pack-out for all seasons
- Patented technology ensures consistent cargo temperatures.
- Meets ISTA (International Safe Transit Association) 3A transit criteria.
- Components are reusable, reducing environmental impact.
- Removable cartridges are pre-conditioned in a freezer.
- Wheels and handles enable easy methods of transport.
- Simple care and cleaning



Cool Cube™ 50

- Lightweight (10 lbs.)
- Easy to carry (includes handle and adjustable shoulder strap).
- Efficient volume: 6" x 6" x 6" (holds 50 multi-dose vials)
- Soft-sided, double zipper case.
- Exterior dimensions: 10.5" x 10.5" x 10.5"
- Pre-condition in 7.25" x 7.25" x 5" of freezer space.

VT-50 (Factory-validated to hold 2° to 8°C for 84+ hours)

VT-50V (Factory-validated to hold -50° to -15°C for 60+ hours)

Cool Cube™ 420

- Large volume: 12" x 12" x 12" (holds 420 multi-dose vials)
- Lockable, airtight case.
- Exterior dimensions: 20.5" x 20.5" x 20"
- Pre-condition in 13.25" x 13.25" x 5.25" of freezer space.
- VT-420 (Factory-validated to hold 2° to 8°C for 101+ hours)
- VT-420V (Factory-validated to hold -50° to -15°C for 113+ hours)
- Easy to pull (includes a collapsible tow handle and four removable swivel wheels).



Response Friendly Medical Products by Vericor

Complete, Organized and Appropriate Quantities for Effective Patient Care

What is Response Friendly?

In essence, Response-Friendly™ is what medical personnel desire in responding to an emergency. When a nurse needs a medium size nitrile glove—it's waiting to be put on. When a surgeon calls for a suture—it's at their fingertips. Everything to handle a medical emergency is where responding personnel expect it to be. It's complete. It's intuitive. It's exclusive. There's only one thing that can replace incident chaos with precise organization and efficiency—a Response-Friendly™ System.

The Response-Friendly™ philosophy is the driving force behind the design and production of VeriCor products. Systems are designed using VeriCor's exclusive software that determines the right quantities for the incident and patient population. Then Flow-Integration™ determines product priority and location. Lastly, equipment and supplies are meticulously organized, comprehensively labeled and intuitively packed to create a Response-Friendly™ System. Only then will a preparedness system perform the way logisticians and health care providers (ultimately patients) need it to—a striking difference from a preparedness stockpile cluttering up a warehouse.

Response-Friendly™ Systems provide:

- Organized storage – ready to satisfy the response time of a preparedness plan.
- Correct equipment – appropriate quantities and compatible components for effective patient care.
- Comprehensive labeling – efficient storage, maintenance, rotation and deployment.
- Product recall monitoring – ensure product integrity.
- A SmartBook™ – a complete system guide; know what you have and where you have it



7 Industrial Way, Salem, NH 03079-2817

evolvetechnology.com

1-888-315-9007