



Emergency Veterinary Equipment

Simple to Use, Easy to Decon, Does Not Sustain Bacterial Growth

All Emergency Veterinary Equipment is constructed of medical grade polyvinyl chloride resin structural tubing and is easily decontaminated. The units will not chip, dent, rust, crack, peel, fade, stain or sustain bacterial growth. These materials allow the units to be very light weight, yet handle a 250 lbs. capacity! The slings are made from durable vinyl coated polyester mesh. The sling and belts can be cleaned with any disinfectant cleaner and can also be removed easily and laundered for deep cleaning.

Stretcher

Our Veterinarian Stretcher can be used by itself or can be combined with one of our Transports.



Elevated Beds

Several sizes and color options available. All beds come with rubber leg covers



Folding Stretcher Transport

Our Veterinarian Transport combined with our Stretcher creates the ultimate modular gurney system. The folding Stretcher Transport is both FOLDABLE & HEIGHT ADJUSTABLE!

Specifications

Folded Dimensions:
5.5"x 26.75"x 63.5"

Height:
28.75" to 35.5"

Width: 26"

Length: 39.5"

Weight: 29 lbs

Weight Capacity:
250 lbs



evolve
emergency
PREPAREDNESS

7 Industrial Way
Salem, NH 03079-2817
1-888-315-9007
evolvetechnology.com



Emergency Veterinary Equipment

Cage Inserts

Custom made to your exact cage dimensions



Quick Adjust Stretcher Transport

Veterinary Transport combined with Stretcher creates the ultimate modular gurney system. The Quick Adjust Stretcher Transport is height adjustable even while the patient occupying it.

Specifications

Weight: 30 lbs.

Width: 26.5"

Height : 36"

Length: 48.5"

Weight Capacity: 250 lbs



Standard Gurney

Standard Gurney is an affordable, easy way to transport injured animals from the family vehicle into the clinic.



Specifications

Gurney Height: 29.5" -39.25"

Width: 26"

Length: 46"

Weight: 38 lbs.

Weight Capacity: 250 lbs.



evolve
emergency
PREPAREDNESS

7 Industrial Way
Salem, NH 03079-2817
1-888-315-9007
evolvetechnology.com